

S-Lock™
Handles &
Blades


Thompson Retractor



- PRECISE BLADE PLACEMENT
- SWIVEL ALIGNS ON TISSUE
- LOCK WHEN SWIVEL IS NOT DESIRED

S-Lock™ Handles

Upgrade your Thompson Retractor frame with S-Lock handles. S-Lock handles provide the same benefits that standard interchangeable handles do but have a lower profile body and contain S-Lock technology, compatible only with the S-Lock blades on the next page.



Handles All S-Lock handle numbers begin with "SL" and require S-Lock blades to use.



S-LOCK CAM II ANGLING HANDLE

Clip-on 8" Angling SL42126WA



S-LOCK 8" HANDLES

8" Angling	SL45003AS
8" Non-Angling	SL45003L
Universal Joint	47204



S-LOCK CAM II HANDLE

Clip-on 8" (non-angling) SL42126WH



S-LOCK BOOKWALTER ADAPTER HANDLES

Angling	SL45005AL
Non-Angling	SL45005L



S-LOCK OMNI ADAPTER HANDLES


Angling	SL45004AL
Non-Angling	SL45004L



S-LOCK SYNFRAME ADAPTER HANDLES

Angling	SL45014AL
Non-Angling	SL45014L


[Not available in Micro-Adjustable handle styles]



Upgrade to Uncompromised Exposure with S-Lock™

New S-Lock Retractor Handles and Blades allow for precise retractor positioning without unwanted movement or rotation to ensure stable retraction in any plane. The swivel locking feature is simple and easy to use, creating a safer operation with reduced time spent repositioning blades.

How it works




1 COMBINE S-LOCK BLADE AND HANDLE

S-Lock handles have serrated faces that match up with serrations on S-Lock blades for precise blade placement.



2 SWIVEL ALIGNS ON TISSUE

Insert S-Lock blade nipple into S-Lock handle by pushing button. Release button with handle on the "swivel" groove nipple setting and swivel to desired position.



3 LOCK WHEN SWIVEL IS NOT DESIRED

Lock S-Lock blade in position on the "lock" groove nipple setting by pushing the button again and pressing straight down over the handle to engage the serrations.

S-Lock™ Blades

Upgrade blades to S-Lock and combine with S-Lock handles to eliminate unwanted blade movement.

All S-Lock blade numbers begin with "SL" and require S-Lock handles to use.



MALLEABLE BLADES

25mm x 110 mm (1" x 4 5/16")	SL45656A
13mm x 90mm (1/2" x 3 1/2")	SL45657A
13mm x 64mm (1/2" x 2 1/2")	SL46115
13mm x 152mm (1/2" x 6")	SL46116
13mm x 203mm (1/2" x 8")	SL46117
25mm x 203mm (1" x 8")	SL46118
25mm x 104mm (1" x 4")	SL46118P
25mm x 254mm (1" x 10")	SL46118A
38mm x 254mm (1 1/2" x 10")	SL46118B
25mm x 127mm (1" x 5")	SL46119
51mm x 102mm (2" x 4")	SL46122P
51mm x 203mm (2" x 8")	SL46122
51mm x 254mm (2" x 10")	SL46122A
38mm x 203mm (1 1/2" x 8")	SL46123
64mm x 254mm (2 1/2" x 10")	SL46124



MALLEABLE FINGER BLADE

6"	SL46191
----	---------



MALLEABLE SLOTTED BLADE

51mm x 203mm (2" x 8")	SL46122V
------------------------	----------



MALLEABLE RENAL VEIN BLADE

25mm x 203mm (1" x 8")	SL46118R
------------------------	----------



MALLEABLE CONCAVE RENAL VEIN BLADES

25mm x 140mm (1" x 5 1/2")	SL46119B
25mm x 191mm (1" x 7 1/2")	SL46119C



HARRINGTON BLADES

38mm x 127mm (1 1/2" x 5")	SL46162P
64mm x 102mm (2 1/2" x 4")	SL46165
64mm x 152mm (2 1/2" x 6")	SL46160
64mm x 203mm (2 1/2" x 8")	SL46163
64mm x 254mm (2 1/2" x 10")	SL46164
114mm x 254mm (4 1/2" x 10")	SL46161



BALFOUR BLADES

30mm x 30mm (1 1/5" x 1 1/5")	SL46142
51mm x 51mm (2" x 2")	SL46142P
76mm x 51mm (3" x 2")	SL46140P



KELLY BLADES

38mm x 51mm (1 1/2" x 2")	SL46127
51mm x 64mm (2" x 2 1/2")	SL46128



DEAVER BLADES (NON MALLEABLE)

25mm x 127mm (1" x 5")	SL46170
38mm x 127mm (1 1/2" x 5")	SL46172
51mm x 127mm (2" x 5")	SL46174



SPLANCHNICH BLADES

51mm x 114mm (2" x 4 1/2")	SL46145
51mm x 152mm (2" x 6")	SL46146



RICHARDSON BLADES

51mm x 76mm (2" x 3")	SL46154
51mm x 102mm (2" x 4")	SL46150



MINI-RICHARDSON BLADES

19mm x 30mm (3/4" x 1 1/5")	SL46125
19mm x 51mm (3/4" x 2")	SL46125A
25mm x 32mm (1" x 1 1/4")	SL46126
51mm x 76mm (2" x 3")	SL46154
44mm x 51mm (1 3/4" x 2")	SL46149



FENESTRATED BLADES

44mm x 89mm (1 3/4" x 3 1/2")	SL46143
48mm x 127mm (1 7/8" x 5")	SL46144

Radiolucent S-Lock™ Blades



RADIOLUCENT MALLEABLE BLADES

25mm x 110mm (1" x 4 1/2")	SL45656AT
13mm x 90mm (1/2" x 3 1/2")	SL45657AT
25mm x 203mm (1" x 8")	SL46118ET
51mm x 254mm (2" x 10")	SL46122AET
51mm x 203mm (2" x 8")	SL46122ET



RADIOLUCENT MALLEABLE CONCAVE RENAL VEIN BLADES

25mm x 140mm (1" x 5 1/2")	SL46119BET
25mm x 191mm (1" x 7 1/2")	SL46119CET



RADIOLUCENT RIGID RENAL VEIN BLADES

25mm x 140mm (1" x 5 1/2")	SL46119RET
25mm x 191mm (1" x 7 1/2")	SL46120RET
25mm x 254mm (1" x 10")	SL46121RET
25mm x 305mm (1" x 12")	SL46122RET



RADIOLUCENT BRAU ANTERIOR LUMBAR BLADES

Reverse Lip ASA Malleable 25mm x 200mm	SL46192ET
Reverse Lip ASA Malleable 32mm x 200mm	SL46193ET
Reverse Lip ASA Rigid 25mm x 100mm	SL46194RET
Reverse Lip ASA Rigid 25mm x 150mm	SL46195RET
Reverse Lip ASA Rigid 25mm x 200mm	SL46196RET
Reverse Lip ASA Rigid 32mm x 100mm	SL46197RET
Reverse Lip ASA Rigid 32mm x 150mm	SL46198RET
Reverse Lip ASA Rigid 32mm x 200mm	SL46199RET
Reverse Lip ASA Rigid 25mm x 255mm	SL46200RET
Reverse Lip ASA Rigid 32mm x 255mm	SL46201RET
Reverse Lip ASA Tapered 25mm x 150mm	SL46195TET
Reverse Lip ASA Tapered 25mm x 200mm	SL46196TET



RADIOLUCENT SPLANCHNICH BLADES

25mm x 152mm (1" x 6")	SL46148ET
25mm x 203mm (1" x 8")	SL46149ET
51mm x 114mm (2" x 4 1/2")	SL46145ET
51mm x 152mm (2" x 6")	SL46146ET
51mm x 203mm (2" x 8")	SL46147ET



RADIOLUCENT DEAVER BLADES

25mm x 127mm (1" x 5")	SL46170ET
38mm x 127mm (1 1/2" x 5")	SL46172ET
51mm x 127mm (2" x 5")	SL46174ET



RADIOLUCENT HARRINGTON

64mm x 152mm (2 1/2" x 6")	SL46160ET
----------------------------	-----------



RADIOLUCENT ODONTOID BLADE

20mm x 90mm	SL45666AT
-------------	-----------



RADIOLUCENT CONCAVE LONGITUDINAL BLADES

23mm x 30mm	SL45160AT
23mm x 40mm	SL45161AT
23mm x 50mm	SL45162AT
23mm x 60mm	SL45163AT
23mm x 70mm	SL45164AT



HEANY BLADES

25mm x 102mm (1" x 4")	SL46180
25mm x 127mm (1" x 5")	SL46182



HOLTGREWE BLADE

60mm x 152mm (2 3/8" x 6")	SL46137
----------------------------	---------



HARRIS PROSTATE RETRACTOR

Prostate Retractor	SL45753AM
--------------------	-----------



DOUBLE ANGLED BLADES

13mm x 70mm (1/2" x 2 3/4")	SL46663
13mm x 76mm (1/2" x 3")	SL46662
25mm x 64mm (1" x 2 1/2")	SL46661
25mm x 76mm (1" x 3")	SL46660



DOUBLE ANGLED BULB BLADE

25mm x 64mm (1" x 2 1/2")	SL46664
---------------------------	---------



VOLKMAN BLADES

Sharp 25mm x 13mm (1" x 1/2")	SL46210
Blunt 25mm x 13mm (1" x 1/2")	SL46211
Sharp 16mm x 13mm (5/8" x 1/2")	SL46212
Blunt 16mm x 13mm (5/8" x 1/2")	SL46213



MURPHY RAKE BLADES

Sharp - 3 Prong 19mm x 13mm (3/4" x 1/2")	SL46216
Sharp - 4 Prong 25mm x 13mm (1" x 1/2")	SL46214
Blunt - 3 Prong 19mm x 13mm (3/4" x 1/2")	SL46217
Blunt - 4 Prong 25mm x 16mm (1" x 5/8")	SL46215



HILL-FERGUSON RECTAL BLADES

32mm x 76mm (1 1/4" x 3")	SL46185
25mm x 64mm (1" x 2 1/2")	SL46186



CONCAVE LONGITUDINAL BLADES

25mm x 30mm (1" x 1 1/4")	SL45160A
25mm x 50mm (1" x 2")	SL45161A
25mm x 70mm (1" x 2 3/4")	SL45162A



CONCAVE MEDIAL BLADES

13mm x 51mm (1/2" x 2")	SL45212A
13mm x 75mm (1/2" x 3")	SL45214A
13mm x 100mm (1/2" x 4")	SL45216A
19mm x 51mm (3/4" x 2")	SL45222A
19mm x 75mm (3/4" x 3")	SL45224A
19mm x 100mm (3/4" x 4")	SL45226A
25mm x 51mm (1" x 2")	SL45232A
25mm x 75mm (1" x 3")	SL45234A
25mm x 100mm (1" x 4")	SL45236A
32mm x 70mm (1 1/4" x 2 3/4")	SL45153A
38mm x 70mm (1 1/2" x 2 3/4")	SL45154A
38mm x 85mm (1 1/2" x 3 3/8")	SL45155A



CONCAVE LATERAL BLADES

9mm x 75mm (3/8" x 3")	SL45640A
19mm x 55mm (3/4" x 2 1/4")	SL45150A
25mm x 30mm (1" x 1 1/4")	SL45156A
25mm x 35mm (1" x 1 1/2")	SL45157A
25mm x 50mm (1" x 2")	SL45151A
25mm x 60mm (1" x 2 3/8")	SL45158A
25mm x 70mm (1" x 2 3/4")	SL45159A
25mm x 85mm (1" x 3 3/8")	SL45152A



CONCAVE TAPERED BLADES

19mm x 35mm (3/4" x 1 3/8")	SL45643A
19mm x 45mm (3/4" x 1 3/4")	SL45645A
19mm x 50mm (3/4" x 2")	SL45646A
45mm x 50mm (1 3/4" x 2 3/8")	SL45641A
45mm x 65mm (1 3/4" x 2 3/4")	SL45642A



"Y" BLADDER BLADE

"Y" Bladder Blade	SL46139
-------------------	---------



ARMY-NAVY BLADES

13mm x 55mm (1/2" x 2 1/8")	SL45620A
13mm x 65mm (1/2" x 2 5/8")	SL45621A
13mm x 75mm (1/2" x 3 1/16")	SL45622A



MIDLINE BLADES

Right 70mm x 83mm (2 3/4" x 3 1/4")	SL45750A
Left 70mm x 83mm (2 3/4" x 3 1/4")	SL45751A
Left 57mm x 90mm (2 1/4" x 3 1/2")	SL45755A
Right 57mm x 90mm (2 1/4" x 3 1/2")	SL45754A



MEYERDING BLADES

16mm x 51mm (5/8" x 2")	SL45614A
16mm x 75mm (5/8" x 3")	SL45600A
16mm x 100mm (5/8" x 4")	SL45601A
25mm x 60mm (1" x 2 3/8")	SL45610A
25mm x 75mm (1" x 3")	SL45611A
25mm x 100mm (1" x 4")	SL45602A
25mm x 125mm (1" x 5")	SL45603A



RADIOLUCENT CONCAVE BLADES

23mm x 35mm	SL46312AT
23mm x 40mm	SL46314AT
23mm x 45mm	SL46316AT
23mm x 50mm	SL46318AT
23mm x 55mm	SL46320AT
23mm x 60mm	SL46322AT
23mm x 65mm	SL46324AT
23mm x 70mm	SL46326AT
23mm x 75mm	SL46328AT
23mm x 80mm	SL46330AT
23mm x 35mm [with teeth]	SL46312ATT
23mm x 40mm [with teeth]	SL46314ATT
23mm x 45mm [with teeth]	SL46316ATT
23mm x 50mm [with teeth]	SL46318ATT
23mm x 55mm [with teeth]	SL46320ATT
23mm x 60mm [with teeth]	SL46322ATT
23mm x 65mm [with teeth]	SL46324ATT
23mm x 70mm [with teeth]	SL46326ATT
23mm x 75mm [with teeth]	SL46328ATT
23mm x 80mm [with teeth]	SL46330ATT



RADIOLUCENT CONCAVE TAPERED BLADES

23mm x 35mm	SL45643AT
23mm x 40mm	SL45644AT
23mm x 45mm	SL45645AT
23mm x 50mm	SL45646AT
23mm x 55mm	SL45647AT
23mm x 60mm	SL45648AT
23mm x 65mm	SL45649AT
23mm x 70mm	SL45650AT
51mm x 50mm	SL45641AT
51mm x 65mm	SL45642AT
23mm x 35mm [with teeth]	SL45643ATT
23mm x 40mm [with teeth]	SL45644ATT
23mm x 45mm [with teeth]	SL45645ATT
23mm x 50mm [with teeth]	SL45646ATT
23mm x 55mm [with teeth]	SL45647ATT
23mm x 60mm [with teeth]	SL45648ATT
23mm x 65mm [with teeth]	SL45649ATT
23mm x 70mm [with teeth]	SL45650ATT
51mm x 50mm [with teeth]	SL45641ATT
51mm x 65mm [with teeth]	SL45642ATT



RADIOLUCENT MIDLINE BLADES

Right 70mm x 83mm	SL45750AT
Left 70mm x 83mm	SL45751AT



RADIOLUCENT CONCAVE MICRO BLADES

14mm x 35mm	SL46510AT
14mm x 40mm	SL46512AT
14mm x 45mm	SL46514AT
14mm x 50mm	SL46516AT
14mm x 55mm	SL46518AT
14mm x 60mm	SL46520AT
14mm x 35mm [with teeth]	SL46510ATT
14mm x 40mm [with teeth]	SL46512ATT
14mm x 45mm [with teeth]	SL46514ATT
14mm x 50mm [with teeth]	SL46516ATT
14mm x 55mm [with teeth]	SL46518ATT
14mm x 60mm [with teeth]	SL46520ATT

[Other blades can be made compatible with S-Lock handles but will only function in the swivel setting. Please call for more information.]



10170 East Cherry Bend Road
Traverse City, Michigan 49684
phone: 231.922.0177
fax: 231.922.0174

thompsonsurgical.com

**For a FREE TRIAL
call 1.800.227.7543**

Free trial valid for U.S. customers only.
International customers please call for more information.